

Introduction

This pincer grip fits to all M-RMJ, M-RMD, M-RMJG and M-RMDG dual shaft sample transporters. Customized jaws can be installed for a variety of different sample holders. Please contact us for a detailed discussion of your requirements.

The opening/closing of the pincer is done with approximately a quarter turn of the outer magnet assembly of the sample transporter.

Description	Part Code
Pincer with customised jaws for sample transporters	M-PGRMS(CUSTOM)

Specifications:

Materials

Stainless steel 1.4301
CuBe2

Vacuum Range

atm to 1×10^{-11} mbar

Bakeout Temperature *

200°C max.

Jaws

Bespoke, as per customer requirements

